ESTEREOMICROSCOPIOS STEREOMICROSCOPES

Serie / Series 200

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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AUXILAB S.L.

CE DECLARATION OF CONFORMITY

STEREOMICROSCOPES SERIES 200 of AUXILAB, S.L. for the Directive of Machines (89/392/CEE modified) and the regulations adopted for their transposition.

NAME OF THE MANUFACTURER/IMPORTER:
AUXILAB, S.L.

ADDRESS: **Polígono Morea Norte, 8
31191 Beriáin (Navarra)**

WE STATE THAT:
**URA TECHNIC STEREOMICROSCOPES SERIES 200
CODES 50200142, 50200141, 50200162, 50200161, 50200182, 50200181, 50200192,
50200191, 50220182 AND 50200201**

Are designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and the national regulations that transpose them.
 - ◆ Directive 73/23/CEE modified over the electric security.
 - ◆ Directive 89/336/CEE modified over the electromagnetic compatibility.
- And that the following harmonized rules have been applied (or part of them):
- UNE 292-1, UNE 292-2, UNE 292-2/A1, UNE 614-1, UNE 1050, UNE 294, UNE 894-1, UNE 894-2, UNE 60204, UNE 61010-1.

BERIAIN 7th July 2004

Signed by: ALFONSO AINCIBURU SANZ
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

return it for overhaul or just to keep it in a safe place when it is not going to be used for a long time.

Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the stereomicroscope and produce an early ageing of the equipment.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives for cleaning the stereomicroscope.
- ◆ Lenses must never be dismantled by the user. In case there is any dirt on the external parts of the lenses, you should clean it with a smooth piece of cloth, fluff-free, dampened with a bit of xylol or toluene.
- ◆ Blow the dust laid down on the lenses with a plastic bulb or clean it either with a soft brush made of natural hair or any special gauze for lenses.
- ◆ For cleaning the mechanical parts you should use non-corrosive lubricants, being careful not to touch the optical parts.

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6. CAUSES OF A DEFECTIVE IMAGE

- In case of a defective image you should check:
- ◆ The illumination is well made.
 - ◆ The candle power is not too weak; if so you should change the lamp.
 - ◆ The cleanliness of the accessible parts of the optical system, by turning the eyepieces and checking that the specks move; if so you should clean them. *You should never dismantle the head*, but you can clean it with delicacy by blowing with a plastic bulb with brush the available parts. If you have an optional objective and the parasite images turn when turning it, you should clean it with a dry brush to remove the dust while you look at it with a magnifying glass or an inverted eyepiece.

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
Aum ⇔ Aumentos

- ◆ Do not place the stereomicroscope near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During its functioning dangerous materials such as flammable or pathological substances must be out of the safety area.
- ◆ When you are not using the stereomicroscope for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ **It is essential to have the equipment switched off and unplugged from the net** before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the stereomicroscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The stereomicroscope is supplied with a Schuko standard wire. It has to be plugged to an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ Never try to repair the stereomicroscope by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...) or damages in nearby equipments.
- ◆ In the event of breakdown please contact your distributor to overhaul through Ura Technic Technical Assistance Department.
- ◆ If the are lamps blown you should replace them by another original Ura Technic ones of 12V10W, being careful not to touch them with the naked hands. Do not use any lamp with higher power, as this could provoke overheating or any other malfunction.
- ◆ Were it essential to replace the fuse (0,5A), you should change it as follows: unscrew the protective lid (FUSE) that you will find at the rear part of the base, replace the fuse locating the new one at the same position and screw it again once you have finished.
- ◆ Do please always use the plastic cover whenever the stereomicroscope is not being used to avoid dust on the optical parts.
- ◆ Please keep the original packaging so as to use it in transport, whenever you have to

♦ The *horizontal arm* is similar in shape to the vertical rod; the first one is placed on the bosshead and tightly screwed with an Allen key of number 2 with two little screws provided. Then you should couple the first ring with the wire. On its end it has a hole in order to assemble an extension to make it flexible.

♦ This *extension* includes a screw to put in the arm's hole and then screw it with the nut with handle without forgetting to put the washers to get a better fastening. At the end of this extension you should put the head, fastening it with two Allen screws of number 2 provided.

♦ The *illuminator* is placed in a circular hole located on the head that is a black piece with cylindrical shape with different widths.

Important: When you are not using the stereomicroscope for a long period of time, please make sure it is disconnected and unplugged from the net.

Setting up

♦ Plug the stereomicroscope to an intake current of 220V 50Hz $\pm 10\%$ provided with an earth connection.

♦ Turn on the most suitable illumination for your needs, incident or transmitted light. There is available a black & white contrast stage with transmitted light as an optional accessory, if necessary.

♦ Centre the sample into the stage and move the focusing control slowly until you manage to see a sharp image observing through the left eyepiece, which is not provided with dioptric compensation ring. Once you have focused the image you should look through the right eyepiece and, if the image is not neat enough, turn the compensation ring until you manage to have a focused image.

♦ With the tool provided (Allen key) you can adjust the focusing control tension to your liking. Thus, you only have to insert an end into the little holes located on the focusing control, maintain the control motionless and turn the key in order to increase/decrease the friction opposed by the focusing control.

♦ For bigger or thicker samples you can adjust the binocular section's height above the base's stand. To do this, you should loosen the rear part screw and, once you have placed it where you need to, tighten it again.

♦ You should observe with both eyes at the same time and graduate the interpupillary distance until the separate images you see become one stereoscopic vision.

♦ Once you have finished working with the stereomicroscope we highly recommend to unplug it, as the electrical components that the adapter has make a little electrical consumption, which produces an early ageing of the equipment and may cause accidents due to heating, as well as an unnecessary waste of energy, etc.

Security

♦ The stereomicroscope must be used by previously qualified staff that knows how the equipment works thanks to the user manual.

♦ You should put the stereomicroscope in a horizontal plane stable table, having a safety area of at least 30 cm per side.

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Installation

- ◆ Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

PLEASE READ THOROUGHLY THE INSTRUCTIONS BEFORE CONNECTING AND OPERATING WITH THIS EQUIPMENT.

- ◆ Please put the stereomicroscope on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the equipment near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.
- ◆ Retire the paper and rubber bands that protect it from dust.
- ◆ Retire also the objective protective cap (model 220).
- ◆ Avoid inflammable or toxic substances in the working area.
- ◆ The stereomicroscope is supplied with a Schuko standard wire
- ◆ Please insert the wire that feeds the AC electric current in the base of current 220V 50Hz $\pm 10\%$ provided with earth wire and to the other end to the stereomicroscope connector.

Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to persons suffered by the improper use of the electric connection. The tension should be 220V 50Hz $\pm 10\%$.

- ◆ You have received a pair of eyepieces together with the stereomicroscope (normally 10x, according to the model supplied). These eyepieces are secured to the tube by means of a screw in order to avoid their loss. Were you to replace them for an optional pair, you should loosen the screw first, then replace them and finally screw them again.
- ◆ If you are not using the stereomicroscope for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

Changing the objectives

- ◆ To change the objectives in models 214 and 216, first you have to loosen two little screws the equipment has on the head, one on the front and the other on the left side. Then, you move the objective to the right until it knocks out of its position, put the desired objective and finally you tighten the screws again.
- ◆ To change the objectives of models 218, 219 and 220 you only have to turn the turret. In model 220 you have to take off the protective cap of the objective before start working.
- ◆ The objective of model 201 is changed in the same way as models 214 and 216. The only thing that changes on the process is that this model only has a screw on the left side of the head instead of having two. This equipment is supplied without assembling, so you have to follow the instructions given below to complete its installation.

Installation of model 201

- ◆ This model has a *base* where are located both the on/off switch and the illuminator plug. It also has a vertical rod with two *rings*: one of them has an inner furrow where the illuminator wire is, and the other one has a screw to adapt the binocular section's height.
- ◆ You have to put a bosshead on the second ring to hold the horizontal arm; this bosshead is meant to adjust height and turn the binocular section.

- Chart 2, models 218, 219 and 220 -

EYEPIECES									
OBJECTIVES	5X		10X		15X		20X		Working Distance
	Mag.	Ø field	Mag.	Ø field	Mag.	Ø field	Mag.	Ø field	
1X	5X	11 mm	10X	20 mm	15X	15 mm	20X	10 mm	57 mm
2X	10X	11 mm	20X	10 mm	30X	7,5 mm	40X	6,5 mm	57 mm
1X	5X	22 mm	10X	20 mm	15X	15 mm	20X	10 mm	57 mm
3X	15X	8 mm	30X	6,7 mm	45X	5 mm	60X	4,3 mm	57 mm
4X	10X	11 mm	20X	10 mm	30X	7,5 mm	40X	6,5 mm	57 mm
6X	20X	5,6 mm	40X	5 mm	60X	3,7 mm	80X	3,2 mm	57 mm

Mag ⇌ Magnifications

- Chart 3, model 201 -

EYEPIECES									
OBJECTIVES	5X		10X		15X		20X		Working Distance
	Mag.	Ø field	Mag.	Ø field	Mag.	Ø field	Mag.	Ø field	
2X	10X	10 mm	20X	10 mm	30X	7,5 mm	40X	6,5 mm	113 mm

Mag ⇌ Magnifications

4. INSTALLATION / SETTING UP

Preliminary inspection

- ◆ Unwrap the stereomicroscope, carefully taking it out and specially taking care of holding it either by the base or the stand but never by the head nor the focusing controls, to avoid focusing precision gears to suffer. Take off the involving plastic and, without connecting the equipment to the net, make sure that it does not present any damage because of the shipment. In case the stereomicroscope presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.
- ◆ Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.
- ◆ Please check that all the accessories are enclosed with the equipment:
 - Objectives (according to model supplied)
 - Ground glass stage
 - 2 lamps 12v and 10w (model 214)
 - 2 spare lamps
 - 2 fuses 0,5 A.
 - Key to adjust the tension (models 220 and 214).
 - Schuko Standard cable.
 - Plastic cover.
 - User manual.
 - Warranty certificate.

We will not accept any equipment in return period unless it comes in its original wrapping.

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3. TECHNICAL SPECIFICATIONS

Ura Technic stereomicroscope's series 200 is composed by the following models and technical specifications:

Model	Objectives	Eyepiece	Lighting		Head		
			Incid.	Transm.	Straight	Inclined 45°	Rotary 360°
210	2X	WF 10X	Yes	No	Yes	No	No
214/1	2X	WF 10X	Yes	No	Yes	No	No
214/2	2X	WF 10X	Yes	Yes	Yes	No	No
216/1	2X	WF 10X	Yes	No	No	Yes	No
216/2	2X	WF 10X	Yes	Yes	No	Yes	No
218/1	2X, 4X	WF 10X	Yes	No	No	Yes	No
218/2	2X, 4X	WF 10X	Yes	Yes	No	Yes	No
219/1	2X, 4X	WF 10X	Yes	No	No	Yes	Yes
219/2	2X, 4X	WF 10X	Yes	Yes	No	Yes	Yes
220	2X, 4X	WF 10X	Yes	Yes	No	Yes	Yes
201	2X	WF 10X	Yes	No	No	Yes	No

Other features (depending on the model supplied)

- Base: metallic and very stable, it contains the transformer for the illumination.
- Pair of slide clips: on the base to hold the simple.
- Binocular section: adjustable in height on the stand.
- Interpupillary distance: adjustable 30°.
- Diopter correction: ± 5 mm in an eyepiece tube.
- Pair of eyepieces: standard (10x), with fixing system.
- Focusing system: by pinion and rack, with safety stop at the end of the route and tension adjustment.
- Focusing controls: bilateral, with friction control system.
- Stage: ground glass of 95 mm of diameter (according to model supplied).
- Pair of rubber eyecups.
- Knob to select the intensity of light: according to model supplied.
- Light switch: selector for incident light, transmitted light or both simultaneously according to the model supplied.
- Halogen lamp 12 V 10W.

Optical specifications

- Chart 1, models 214 and 216 -

OBJECTIVES	EYEPIECES								Working Distance
	5X		10X		15X		20X		
	Mag.	Ø field	Mag.	Ø field	Mag.	Ø field	Mag.	Ø field	
1X	5X	22 mm	10X	20 mm	15X	15 mm	20X	10 mm	57 mm
2X	10X	12 mm	20X	10 mm	30X	7,5 mm	40X	6,5 mm	80 mm
3X	15X	8 mm	30X	6,7 mm	45X	5 mm	60X	4,3 mm	61 mm
4X	20X	6 mm	40X	5 mm	60X	3,7 mm	80X	3,2 mm	57 mm
6X	30X	4 mm	60X	3,3 mm	90X	2,5 mm	120X	1,8 mm	66 mm

Mag ⇄ Magnifications

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT.....	12
2. DESCRIPTION	12
3. TECHNICAL SPECIFICATIONS	13
4. INSTALLATION / SETTING UP.....	14
5. MAINTENANCE AND CLEANING	17
6. CAUSES OF A DEFECTIVE IMAGE	18
ANNEX I: CE CERTIFICATE.....	19

1. USES OF THE INSTRUMENT

Stereomicroscopes are used to observe enlarged images of little objects, being the image you see through it a stereoscopic one, not inverted and preserving the object's original colour.


As a result of this, they can be used in a wide range of fields going from education to electronic industry, passing through every single activity that needs an object to be seen enlarged but maintaining a stereoscopic (three-dimensional) vision.

We stand out some field where stereomicroscopes can be used:

- Education: primary, secondary, vocational training and university.
- Industry: electronics, clock making, textile, jewellery, painting, food industry, tool setup, casting (metallurgy), graphics, galvanising, etc. In every process of inspection, assembly and repairing of instruments, measuring plates or any other kind of precision piece, as well as in quality control processes.
- Agriculture and stockbreeding, biology, prosthodontists, artistic restoration, philately, mycology, graphology, mineralogy, optics, research, etc.

These equipments provide a great comfort in their use and an excellent quality of the observed image.

2. DESCRIPTION

1. Holding clips
2. Stage
3. Base
4. On/off switch
5. Incident light
6. Focusing control knob
7. Height adjustment knob for binocular section
8. Focusing tension adjustment ring
9. Eyepieces
10. Dioptic adjustment ring
11. Head
12. Objective/s
13. Incident and transmitted light

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Thank you for choosing this equipment. We sincerely wish that you enjoy your stereomicroscope Ura Technic series 200. We highly recommend looking after this equipment according to what is stated in this manual.

Ura Technic develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of *Maintenance*, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ This manual is inseparable from the stereomicroscope Ura Technic series 200, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the stereomicroscope avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the stereomicroscope to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the equipment.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Ura Technic Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0,30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the stereomicroscope warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the stereomicroscope warranty.
- ◆ Neither fuses (0,5A) nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Ura Technic Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!!
IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

INDEX OF LANGUAGES

Spanish	2-10
English	11-19

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958 412 886

629 636 705

<http://www.cromakit.es/>

Calle Tucumán 8 Nave B 18200 Maracena (Granada)