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Directive 2002/96 / EC

This symbol on the nameplate of the equipment or on its packaging indicates that this product may not be treated as household waste and therefore must be disposed of at a collect point of electrical and electronic equipment.

With this action, you. You are helping to better preserve the environment.

For more detailed information about recycling of this equipment, please contact your city office, your household or nearest dealer where you purchased your equipment.

## INTRODUCTION & GUARANTEE

Welcome to the bath for education RS 3.5 user manual. Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

When help needed, you can always contact the service department of manufacturer or your supplier for technical support. Please provide the customer care representative with the following information: Serial number, description of problem, methods and procedures adopted to resolve the problems and your contact information.

These baths are warranted to be free from defects in materials and workmanship under normal use and service for a period of 12 months from the date of invoice. The breakdowns caused by inadequate use, handling errors or any other reason not attributable to a manufacturing defect are excluded from this guarantee.

For claims under the warranty please contact your local dealer. You may also send the instrument direct to our works, enclosing the invoice copy and giving reasons for the claim.

## FEATURES

- Power supply; 230V-50Hz
- Equipment designed for indoor use.
- Control box with metal anodized finish.
- Transparent polycarbonate tank capacity of 3.5L.
- Heating by means of electric resistance shielded in stainless steel AISI 304
- Temperature regulation by means of bulb thermostat from ambient +5°C to 90°C.
- Equipped with safety thermostat with manual reset, in case of over-heating or lack of fluid, according to DIN 12877 class 2
- Stainless steel rack for 18 tubes up to Ø20mm.
- Stainless steel support to prevent direct contact of the steam with the control box

<i>Code</i>	57 000 035 / 50600035
<i>Power</i>	1000 W
<i>Capacity of the tray</i>	3.5 L
<i>Overall dimensions of the tray</i>	160 x 150 x 264 mm (W x H x D)
<i>Dimensions with control box &amp; stand</i>	160 x 225 x 264 mm (W x H x D)
<i>Maximum working temperature</i>	90°C

## INSTALLATION

Equipment must be placed on to flat and robust surface, without any ground level differences.

## OPERATION

Control panel :

1. Fuse.
2. Thermostat security operation pilot light (red).
3. Heating pilot light.
4. Luminous operation switch.
5. Thermostat control.
6. Safety thermostat.

For your own safety, do not switch the bath on if it is not properly earthed.

Before connecting the apparatus to the mains, make sure the voltage that appears on the characteristics plate is the same as the one from the mains.

Fill the bath up with liquid, making sure always that the level of liquid is above the resistor of the apparatus.

Connect the apparatus to the mains.

Place the operation switch on the I position.

Turn the thermostat control up to the desired working temperature. The heating pilot light will light up at the same time.

While heating the liquid, the amber pilot light will light up to show resistor operation. Once the desired operating temperature has been reached the amber light will go off and only come on and off intermittently to maintain the said temperature.

Keep evaporation of the liquid in mind so that the resistor of the bath is never left bare.

The graduation of the scale of the thermostat, although close to the real one, it is always a rough guide. The real temperature is the one that will be indicated by a thermometer placed on the surface of the bath.

## SAFETY RECOMMENDATIONS

Follow the safety recommendations to prevent damage to persons or property.

- Connect the unit only to a grounded mains power socket.
- Operation is permitted with non-flammable liquids only.
- Never operate the unit without bath fluid in the bath.
- Be careful when emptying hot bath fluids. Check the temperature of the bath fluid before draining.
  - Never operate damaged or leaking equipment.
  - Always turn off the unit and disconnect the mains cable from the power source before performing any service or maintenance procedures or before moving the unit.
  - Always empty the bath before moving the unit.
  - Never operate equipment with damaged mains power cables.
- Some parts of the cover bath may become extremely warm during continuous operation. Therefore, exercise particular caution when touching these parts.

The device includes a safety thermostat, on the rear side of the bath.

Its use is to protect the element in the event of over-heating, normally due to a lack of liquid in the bath.

In the event this is activated, the red light will be switched on and the element will be disconnected avoiding the operation of the bath (does not heat up).

In order to start the bath again, you have to wait for the temperature to decrease by around 50%; then press the re-arm switch on said thermostat (small red button).

## CLEANING THE UNIT

Before cleaning the unit, disconnect the power plug from the mains socket.

Sometimes it's happens that if water of the network is used, depositions of lime in the bath tank and in the resistance be produced. You must avoid that the resistance work cover with lime, because the resistance heats less and its useful life diminishes.

To clean of lime the bath tank to fill it of water with a 10% rate of nitric acid, until all the lime be removed. After empty it and rinse it with clean water

If you don't have nitric acid you can use hot water mixed with two glass of wine vinegar.

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958 412 886

629 636 705

<http://www.cromakit.es/>

**Calle Tucumán 8 Nave B 18200 Maracena (Granada)**