**CK-MB Control serum  
IVD**

Store at 2 - 8°C

**PRODUCT CHARACTERISTICS**

CK-Nac / CK-MB CONTROL is a human lyophilised control serum.

Is intended for control of accuracy for use with manual and automated analytical procedures.

**REAGENTS**

Human serum. Biological additives. Bacteriostatics agents. The concentration of the components are lot-specific. See chart below for the values of ingredients.

**PRECAUTIONS**

Components from human origin have been tested and found negative for the presence of HBsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**PREPARATION**

-Reconstitute (→) with 2.0 mL of distilled water.  
-Mix thoroughly, avoiding foam forming.  
-Bring to room temperature for about 30 min. before use. Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

**STORAGE AND STABILITY**

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth. Store tightly capped when not use.

- After reconstitution is stable for:
- At 15°C to 25°C 24 hours
- At 2°C to 8°C 3 days
- At -25°C to -15°C 1 month

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