

▶ VORTEX STIRRER RSLAB-6

- 1| Suitable for touch or continuous operation.
- 2| Adapters securely clicked onto appliance in any position.
- 3| Sturdy cast casing.
- 4| Specially designed silicone feet and steel base for excellent vibration damping.
- 5| Eccentric with oil-less ball bearings.
- 6| Protection class according to DIN EN60529: IP21.
- 7| Provided with a standard top for tubes with diameter <30 mm. (code 90681581)




model RSLAB-6



model RSLAB-6PRO

▶ technical features

Reference	50681001	50681501
Model	RSLAB-6	RSLAB-6PRO
Speed range	500 rpm (fixed)	0-2500 rpm
Shaking movement	Orbital (diameter 4mm)	
Dimensions	127x130x160 mm	
Weight	3,5 Kg	
Protection class	IP21	



STIRRERS

VORTEX

RS Lab

▶ Accessories for RSLAB-6PRO model (Not included)



Reference

90681511

Universal top for tube adapters from code 90681521 to 90681591.



Reference

90681571

Tube holding rod from code 90681521 to 90681561.



Reference

90681521*

Adapter for 48 PCR microtubes 6mm.



Reference

90681531

Adapter for 15 round bottom tubes 10 mm.




Reference

90681541

Adapter for 12 round bottom tubes 12 mm.

*For use in continuous mode.



STIRRERS

VORTEX

RS *Lab*



Reference

90681551

Adapter for 8 round bottom tubes 16 mm.



Reference

90681561*

Adapter for 8 round bottom tubes 20 mm.



Reference

90681591

Antiskid platform for tubes max. diameter 99 mm.




Reference

90681582

Rubber feet.

*For use in continuous mode.



Distribuido por



958 412 886



629 636 705



<http://www.cromakit.es/>

Calle Tucumán 8 Nave B 18200 Maracena (Granada)

